VISKO, VISKODU, VISKOTRYS, KUBA, MEDAS

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EN This Instruction for use **MUST** be given to the user of the product.
BEFORE using this product, this Instruction for use **MUST** be read and saved for future reference.

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SYMBOLS

In this Instruction for Use warnings are indicated by symbols.
The warning symbols are accompanied by a heading that indicates the severity of the danger or basic information.

Gives basic information, recommendations and information for efficient, trouble-free use.

WARNING / CAUTION - to indicate that caution is necessary when operating the device. To indicate that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.

ATTENTION! This sign is used for very important parts concerning your safety. Please look at the operating instructions.

SAFETY INFORMATION

WARNING
For indoor use only!

WARNING
Keep away from fire!

CAUTION
It's recommended to place the mattress on the solid base, which would be air permeability.

CAUTION
Mattress has an elastic component, so there can appear permissible 1-2 cm deflection of the size. The properties of the foam in the places which have a load might be changed. This property is not considered as a spoilage.

CAUTION
Do not use a knife to open the package, you can damage the mattress.

GENERAL

INTRODUCTION

This Instruction for use contains important information about the handling of the product. To ensure safety when using the product, read the Instruction for use carefully and follow the safety instructions.

Medical mattress VISKO / VISKODU / VISKOTRYS / KUBA / MEDAS is non-allergenic, breathable, washable and autoclavable as well as meets medical standards. The package includes a medical cover that elastically covers the mattress on all sides and can be removed by zippers. The fabric is liquid-resistant, hygienic, breathable material. The unique properties of the fabric do not allow bacteria to multiply, the cover is fireproof and has a long service life.

This instruction for use is available in multiple languages to comply with MDR 2017/745. Please ensure that the correct language version is provided to the user in their respective Member State.

i Note that there may be sections in this document, which are not relevant to your product, since this document applies to all available models. If not otherwise stated, each section in this document refers to all models of the product.

PREFACE

Thank you for the confidence shown in Kid-Man. Every product is factory checked and left us in faultless condition.

We wish you much success and we hope UAB Kasko Group | Kid-Man® product will make a positive impact to your daily living.

WARRANTY

We issue a warranty of 24 months from date of purchase. Any unauthorized modifications carried out without approval of the manufacturer will make the liability null and void. This includes exceeding the maximum user loads and using product in improper way. Warranty is only valid with fully filled warranty card. UAB „Kasko Group“ after getting claimed item from purchaser within period of 4 weeks decides whether issue is treated under warranty. If issue is not covered under warranty and customer agrees, UAB „Kasko Group“, if possible, fixes issue. All costs, including costs of transportation, should be covered by customer. Warranty repair does not extend warranty.

LIMITATION OF LIABILITY

Manufacturer accepts no liability for damage arising from:

- Non-compliance with the Instruction for use;
- Incorrect use;
- Natural wear and tear;
- Incorrect assembly or set-up by the purchaser or a third-party technical modifications;
- The usage of unapproved 3rd party accessories instead of approved Manufacturer's accessories;
- Unauthorized modifications and/or use of unsuitable spare part.

STANDARDS AND REGULATIONS

Medical mattress complies the requirements of:

a) REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 5 April 2017 on medical devices, amending Directive 2001/83/EC, Regulation (EC) No 178/2002 and Regulation (EC) No 1223/2009 and repealing Council Directives 90/385/EEC and 93/42/EEC.

ISO 14971:2019 Medical devices — Application of risk management to medical devices ISO 14971:2019

ISO 10993-1:2018 Biological evaluation of medical devices — Part 1: Evaluation and testing within a risk management process ISO 10993-1:2018.

ISO 15223-1:2021 Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements ISO 15223-1:2021.

SERVICE LIFE

The expected service life is 5 (five) years, presuming that the product is used daily and in accordance with safety instructions, maintenance instructions and intended use, stated in this Instruction for Use. After then, medical mattress has to be recycled. Do not throw your medical mattress away with common waste (see chapter “DISPOSAL”).

INFORMATION ABOUT MANUFACTURER

Manufacturer: UAB “Kasko Group” (registered trade mark Kid-Man®)

Address: A. Vivulskio 12D, Vilnius, LT-03221, Lithuania

Phone: +370 5 265 0000

E-mail: info@kid-man.com

Website: www.kid-man.com

INFORMATION ABOUT PRODUCT

Product: Medical mattress

Model: VISKO, VISKODU, VISKOTRYS, KUBA, MEDAS

REF: KUBA, MED200901025, MED200901035, MED200901225, MED200901235, MED200901225C, MED200901425, VISKO-2, VISKO-2/20010015, VISKO-2/20012015, VISKO-3, VISKO-3/20012015, VISKO, VISKO/2009012

INTENDED PURPOSE AND INDICATIONS

These comfortable medical mattresses are designed for patients with different levels of bedsores risk relief needs. See the table (Table 1) below to identify the intended purpose of your medical mattress according to your model:

Table 1

Medical mattress model	MEDAS	KUBA	VISKO	VISKODU	VISKOTRYS
Bedsores risk relief level	Very low	Low	Medium	High	Very high

The mattress is not intended for use in environments with extreme temperature variations outside the specified operating range (0°C to 50°C) or in applications beyond its stated purpose.

INTENDED USERS

MEDICAL MATTRESS is intended to be used by the patients with a risk of bedsores. The medical mattress model is adjusted to the specific needs of the patient (see Table 1).

CONTRA-INDICATIONS

There are no contraindications known, when using the medical mattress as intended.

RESIDUAL RISKS

Medical mattresses designed to prevent and help to treat bedsores (pressure ulcers) are essential in patient care, especially for individuals with limited mobility. However, even with advanced mattress technology, certain residual risks remain:

1. Incomplete Pressure Redistribution:

While specialized mattresses aim to alleviate pressure on vulnerable areas, they may not entirely eliminate it. Prolonged pressure, even at reduced levels, can still lead to tissue damage in susceptible patients.

2. Shear and Friction Forces:

Mattresses may not fully prevent shear (the sliding of skin layers) and friction, which are significant contributors to bed sore development. Patients who frequently adjust their position or are repositioned without proper technique remain at risk.

3. Moisture Accumulation:

Some mattress materials can retain moisture, creating a damp environment that softens the skin and increases susceptibility to ulcers. This is particularly concerning for patients with incontinence or excessive sweating.

4. Heat Retention:

Certain mattress types, such as foam mattresses, may retain body heat, leading to increased skin temperature and perspiration, which can exacerbate skin breakdown.

5. Inadequate Individualization:

A mattress that is not appropriately selected for a patient's specific needs—considering factors like body weight, medical condition, and mobility level—may not provide optimal protection, leaving certain areas at risk.

6. Damaged or Worn Covers:

Over time, mattress covers can become damaged or worn, compromising their protective barrier and increasing the risk of infection. Regular inspection and maintenance are crucial to mitigate this risk.

The identified residual risks have been evaluated as acceptable within the risk-benefit analysis conducted in accordance with ISO 14971. Comprehensive risk management measures have been implemented to minimize potential harm.

SIDE EFFECTS

Medical mattresses designed for patients at risk of bedsores (pressure ulcers) are essential in healthcare settings. However, their use can sometimes be associated with certain side effects:

1. Skin Maceration:

Some mattresses may retain moisture, leading to skin maceration, where the skin becomes overly soft and prone to breakdown. This is particularly concerning for patients with incontinence or excessive sweating.

2. Heat Retention:

Certain mattress materials can retain body heat, potentially causing increased skin temperature and perspiration, which may exacerbate skin breakdown.


3. Allergic Reactions:

Some patients may experience allergic reactions to materials used in mattresses, such as latex or specific synthetic compounds, leading to skin irritation or dermatitis.

PRODUCT LABEL

The identification label is attached to the product (see picture).

1. Product
2. Model
3. Maximum load
4. Manufacturer / contact details
5. Date of manufacture
6. UDI / SN code
7. CE mark
8. Consult instructions for use

CONTENT OF PACKAGE


The medical mattress is delivered rolled up in a package. Once unrolled, please keep mattress flat without load for 24 hours before use. (refer to page 10). **Warning!** Using mattress immediately after unpacking it may cause harm of foam insert. After unpacking, please keep the packaging if possible. You may need it later for storage or for the return of the mattress. Once parcel is received, please check whether it is not damaged during transportation.

i If you see signs of damage, please make note to transport company on the same time.

i Please check content of package and if any of below parts are missing, get in touch with your local dealer.

Content of package:

1. 1 x mattress
2. 1 x mattress protective cover
3. Instruction for use


COMPONENTS AND FUNCTIONS

VISKO

Overview


1. Waterproof cover with U shape zipper
2. Memory foam **60 kg/m³**
3. Standard foam **38 kg/m³**

VIŠKODU

Overview


1. Waterproof cover with U shape zipper
2. Memory foam **60 kg/m³**
3. High resilience foam **42 kg/m³**

VIŠKOTRYŠ

Overview


1. Waterproof cover with U shape zipper
2. Memory foam **60 kg/m³**
3. High resilience foam **37 kg/m³**
4. High resilience foam **42 kg/m³**

KUBA

Overview


1. Waterproof cover with U shape zipper
2. High resilience foam **35 kg/m³**

MEDAS

Overview

1. Waterproof cover with U shape zipper
2. Standard foam **25 or 35 kg/m³**

i Complexations of your medical mattress can differ a little from a medical mattress described in this manual due to continuous development of the product.


COMPONENTS AND FUNCTIONS

- Read the instruction for use (Picture 1 (1)).
- Open the package using scissors (Picture 1 (2)).
- Unroll the medical mattress (Picture 1 (3)).
- Let the mattress fully expand for 24 hours (Picture 1 (4)).

CAUTION

Do not use a knife to open the package, you can damage the mattress.

Picture 1

OPERATION

PLACING THE MEDICAL MATTRESS

It's recommended to place the mattress on the solid base, which would be air permeability.

CAUTION

Mattress has an elastic component, so there can appear permissible 1-2 cm deflection of the size. The properties of the foam in the places which have a load might be changed. This property is not considered as a spoilage.

MAINTENANCE

To ensure safe and reliable operation, carry out the following checks and maintenance regularly or have it carried out by another person. Below You will find recommended maintenance's schedule.

	Before use	Once a week	4-6 time per year	Once a year
Visual check	X			
Ventilate the mattress			X	
Clean the mattress cover		X		
Have medical mattress checked by an authorized dealer				X

WARNING

Repair or adjustments of medical mattress should be performed by authorized service person while person is not using the product.

CLEANING AND DISINFECTION

Mattress cover cleaning can be performed by disinfection in autoclave at +120 °C, or washing up at +95 °C, and can be dried in a dryer at up to +70 °C.

Foam can be also disinfected in a disinfection chamber for 15-20 min at 105 °C temperature.

Mattress can be ventilated regularly (recommended 4-6 times a year).

Cleaning and disinfection can be performed by following the instructions given by the manufacturer.

STORAGE

The medical mattress should be stored in a dry environment, away from direct sunlight. The foam mattress is resistant to temperature from 0°C to +50°C.

RE-USE

The medical mattress is suitable for re-use. The following actions should be taken:

Cleaning and disinfection (see chapter „CLEANING AND DISINFECTION“); Inspection (see chapter „MAINTENANCE“).

DISPOSAL

Be environmentally aware and dispose of your medical mattress properly. Disposal is subject to national and local regulations.

To have medical mattress properly recycled, please contact your local community, disposal facility or our service department.

SPECIFICATIONS

Model	VISKO	VISKODU	VISKOTRYS	KUBA	MEDAS 25 kg/m ³	MEDAS 35 kg/m ³
Maximum safe load (kg)	140	200	240	110	120	120
Foam insert weight (kg)	12	13	14	7	4,5 / 5,4 / 6,3	6,3 / 7,6
Cover weight (kg)	0,5	0,5	0,5	0,5	0,5	0,5
Overall height (cm)	14	15	15	11	10 / 12 / 14	10 / 12
Overall width (cm)	90	90	90	90	90	90
Overall length (cm)	200	200	200	200	200	200
Zipper shape of waterproof cover	“U” shape	“U” shape	“U” shape	“U” shape	“U” shape	“U” shape

INCIDENTS

Any serious incident that has occurred in relation to the device should be reported to UAB Kasko Group | Kid-Man® and the competent authority of the Member State in which the user and/or patient is established.

CONTACTS:

Manufacturer: UAB “Kasko Group” (registered trade mark Kid-Man®)

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