

KYLA

Lifting mat

Key features & benefits:

- KYLA lifting mat is designed for positioning, transferring a lying patient from one surface to another, and lifting or turning the patient in bed.
- Available with 4 or 6 reinforced handles to support controlled multi-person handling and allow selection according to transfer needs.
- Reinforced handle layout helps caregivers maintain a secure grip, improving control, weight distribution and safety during transfers.
- The flexible, one-piece plastic board has foam padding on one side, providing stable support and improved user comfort during movement.
- Lightweight construction makes the lifting mat easy to position, handle, store and use in daily care routines.
- Two size options allow selection according to user height and required transfer support.
- Durable PU-coated rayon fabric supports hygiene and regular care.
The 64% rayon / 36% polyurethane composition is specially selected for higher washing temperature requirements, making this an improved version for more demanding care environments.
- Material properties support professional use, including good tensile strength, tear strength, peel adhesion and hydrostatic pressure resistance for repeated use. (Before washing the red cover, open the Velcro closure and remove the padded flexible plastic insert from inside the mat.)
- The double-layer removable cover is washable at 60 °C, supporting hygienic and long-term use.



06-4703



06-4704

Product attributes:



Washable
cover 60 °C



Light-
weight

Colour range:



Main information:

REF	Size	Weight, kg	Outside material	Inside material
06-4703	110 x 53	1.2	64% Rayon 36% Polyurethane	Flexible plastic board with foam padding
06-4704	150 x 53	1.6		

Technical data:

Overall thickness, mm	General width, cm	Padded plastic board width, cm	General length, cm	Max Load, kg
5	53	31.5	110 / 150	150 / 180

Visuals:



Maintenance:

It is recommended to store the lifting mat in a dry place, protected from direct sunlight. Maintain it regularly according to the instructions on the product label. Do not use abrasive cleaners. The flexible plastic base can be washed and disinfected.

Warranty:

The device has a 24-month warranty from the date of purchase.

Packaging / logistic information:

Packaging / logistic information:

PER BOX

PER PALLET

REF	BOX DIMENSIONS	PER BOX	PER PALLET	EAN
06-4703	L1100 x W340 x H30 mm	1 pc	110 boxes	4 779031 631816
06-4704	L1500 x W340 x H30 mm	1 pc	100 boxes	4 779031 631809

Storage conditions:



Manufacturer information:

 **UAB Kasko Group | Kid-Man®**
A.Vivulskio str. 12D, 03221, Vilnius, Lithuania
www.kid-man.com | info@kid-man.com
Tel nr: +370 5 265 0000

Importator / Distributor in Romania : Kid-Man Aid SRL
Strada Theodor Aman Nr 2-2A, Braila Tel nr: +40 753 660066

