


# Kid-Man<sup>®</sup> Examination gloves

## TECHNICAL DATA SHEET





## Key features of NITRILE GLOVES

- Thin powder-free nitrile glove, ambidextrous and with rolled rim for extra cuff strength
- Good viral barrier properties, protects against a variety of chemicals
- Excellent tactile sensitivity, enhanced wet and dry grip with microtextured fingertips
- Extremely tear-resistant
- Excellent fit, easy to don and remove
- Free of latex, powder and proteins, good alternative to those suffering from type I allergies
- The gloves are suitable for all food contact



MEDICAL	2017/745	CLASS I	EN 455	AQL 1.5
PPE	2016/425	CATEGORY III	EN 374	AQL 1.5



NIGL100AF (XS-XL)

### Accelerator-free For allergic skin

- Sulphur-free, chemical accelerator-free, latex-free and powder-free
- Skin friendly and reduces the risk of allergies
- Beaded cuff ensures easy donning and prevents rolling down
- Comfortable and reliable hand protection

#### Application areas

- Chemical industry
- Workshop activities
- Painting, varnishing and polishing work
- Electronics and computer industry
- Laboratory and research
- Cleaning and maintenance
- Product protection and printing industry



NIGL100290 (XS-XL)

### 290 mm length More safety & comfort

- Longer and stronger
- Sulphur-free, latex-free and powder-free
- Skin friendly and reduces the risk of allergies
- Beaded cuff ensures easy donning and prevents rolling down
- Comfortable and reliable hand protection

#### Application areas

- Chemical industry
- Workshop activities
- Painting, varnishing and polishing work
- Electronics and computer industry
- Laboratory and research
- Cleaning and maintenance
- Product protection and printing industry



NIGL50GRIP (M-XXL)

### Diamond texture For optimal grip

- Extremely tear-resistant
- Silicone-free manufacturing process
- Reliable hand protection

#### Application areas

- Heavy-duty work, vehicle repair
- Chemical industry
- Workshop activities
- Painting, varnishing and polishing work
- Electronics and computer industry
- Laboratory and research
- Cleaning and maintenance
- Product protection and printing industry



NIGL100BDG (XS-XL)

### Biodegradable Eco-friendly choice

- Biodegradable gloves
- Biodegradation within 270 days
- Contains innovative biodegradable technology
- Comfortable and reliable hand protection

#### Application areas

- Green and sustainable companies
- Chemical industry
- Workshop activities
- Painting, varnishing and polishing work
- Electronics and computer industry
- Laboratory and research
- Cleaning and maintenance
- Product protection and printing industry



## Key features of VINYL & LATEX gloves



LTGL100 (XS-XL)

### Latex

- Made of natural rubber (NR) latex, which offers good resistance against many acids and alkalis
- The best fit and feel, comfortable to wear
- A disadvantage of natural latex is that its natural proteins can cause or trigger latex allergies.
- Protein content does not exceed 50 µg/g.

#### Application areas

- Healthcare, nursing and care homes
- Dental clinics
- Beauty and tattoo studios
- Food industry and catering
- Laboratory and research
- Cleaning and maintenance
- Household and janitorial work



VIGL100 (S-XL)

### Vinyl

- Thin powder-free glove, ambidextrous and with rolled rim for extra cuff strength
- Good viral barrier properties, protects against a variety of chemicals
- Easy to don and remove
- Free of latex, powder and proteins, good alternative to those suffering from type I allergies
- The gloves are suitable for all food contact.

#### Application areas

- Food industry and catering (popular among restaurant, takeaway and sushi bar kitchens)
- Beauty and hair salons
- Cleaning and maintenance
- Retail and cashiers
- Light workshop activities
- Household and janitorial work
- Product protection and packing

## Technical details about Kid-Man® gloves assortment

### Size range

What sizes of gloves can you get?






### Colour range




What colours of gloves can you get?



### Which glove material should I choose?

Individuals may perceive the material properties of nitrile, vinyl and latex differently. The following comparison provides a general overview and summarises the main material characteristics to help you choose the right type of glove.

Does it match what you are looking for?  The best choice  Recommended  Not recommended

Material	NITRILE								VINYL	LATEX	
	REF / Product name	NIGL100 NIGL100BL NIGL100RO	NIGL100BK	NIGL100290	NIGL100AF	NIGL100BDG	NIGL100OAT	NIGL100EKO	NIGL50GRIP	VIGL100	LTGL100
Comfort		✓	✓	✓	✓	✓	✓	✓	✓	✗	
Elasticity		✓	✓	✓	✓	✓	✓	✓	✓	✗	✓
Grip / Tactility		✓	✓	✓	✓	✓	✓	✓		✗	✓
Tear resistance		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Elongation		600 %	600 %	600 %	600 %	600 %	600 %	600 %	600 %	300 %	800 %
Puncture resistance		✓	✓	✓	✓	✓	✓	✓		✓	✓
Chemical resistance		✓	✓	✓	✓	✓	✓	✓	✓	✗	✓
Latex-free		✓	✓	✓	✓	✓	✓	✓	✓	✓	✗
Accelerator-free		✗	✗	✗		✗	✗	✗	✗	✓	✗
Plasticisers		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Heavy duty / extra thick		✗	✗	✗	✗	✗	✗	✗		✓	✗
Ecological		✗	✗	✗	✗		✗	✗	✗	✗	✗
Size range		XS-XL	XS-XXL	XS-XL	XS-XL	XS-XL	XS-XL	XS-XXL	M-XXL	S-XL	XS-XL
Colour		Violet Blue Pink	Black	Blue	Blue	D. Blue	Blue	Blue	Black	Clear	White

### What are the logistics of transporting the gloves?

Sales unit: transport carton.



REF / MODEL NAME	GLOVES PER PACKAGE (PCS.)	PACKAGES PER TRANSPORT CARTON (PCS.)	PACKAGES ON THE PALLET (PCS.)
NIGL100, NIGL100BK, NIGL100BL, NIGL100RO	100 (XXL size 90 pcs.)	10	1080
NIGL100AF, VIGL, LTGL	100	10	800
NIGL100290	100	10	630
NIGL100BDG	100	10	840
NIGL50GRIP	50	10	800

## MAIN FEATURES

- Thin powder-free nitrile glove, ambidextrous and with rolled rim for extra cuff strength
- Good viral barrier properties, protects against variety of chemicals
- Excellent tactile sensitivity, enhanced wet and dry grip with microtextured fingertips
- Excellent fit, easy to don and remove
- Free of latex, powder and proteins, good alternative to those suffering from type I allergies
- The gloves are suitable for all food contact



Color:  
Violet

### Recommendations for use:



Basic Procedures and wards



Cleaning



Food Handling

MEDICAL	2017/745	CLASS I	EN 455	AQL 1.5
PPE	2016/425	CATEGORY III	EN 374	AQL 1.5

## MAIN INFORMATION

SIZE	REF	PALM WIDTH	OVERALL LENGTH	THICKNESS		
				PALM	CUFF	FINGER
XS	NIGLXS100	75 ± 5 mm	≥240 mm	Min 0.08 mm	Min 0.07 mm	Min 0.10 mm
S	NIGLS100	85 ± 5 mm				
M	NIGLM100	95 ± 5 mm				
L	NIGLL100	110 ± 5 mm				
XL	NIGLXL100	120 ± 5 mm				

MATERIAL	EXTERIOR	INSIDE TREATMENT	FORCE AT BREAK		SHELF TIME	NATURAL LATEX FREE	ACCELERATOR FREE
			Before ageing	After ageing			
Nitrile	Fingertip textured	Powder Free	≥6N		5 years	Yes	No

## LABELLING

- Each unit displays the following information: product reference, batch number and expiry date.
- Symbols used in compliance with standard EN ISO 15223: Symbols to be used with medical device labels, labelling and information to be supplied.

**CE 2777**

- Notified Body responsible for certification and ongoing conformity: SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland



**EN ISO 374-5:2016**  
Protects against viruses

- Resistance to bacteria and fungi - pass. Resistance to virus - pass.\*The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

## PACKAGING / LOGISTIC INFORMATION

AVAILABLE SIZES	PCS/DISPENSER	PCS/TRANSPORT CARTON	DISPENSER	TRANSPORT CARTON
XS	100 pcs.	1000 pcs.		
S	100 pcs.	1000 pcs.		
M	100 pcs.	1000 pcs.		
L	100 pcs.	1000 pcs.		
XL	100 pcs.	1000 pcs.		



50 PCS PER BOX



10 PCS PER BOX



1080 BOXES PER PALETTE

**Sales Unit:** Transport Carton

**Storage conditions:**




**UAB Kasko Group | Kid-Man®**

A.Vivulskio str. 12D, 03221, Vilnius, Lithuania  
www.kid-man.com | info@kid-man.com  
Tel. +370 5 265 0000

Importator / Distributor RO : Kid-Man Aid SRL  
Strada Theodor Aman Nr 2-2A, Braila Tel: +40 753 660066

## MAIN FEATURES

- Thin powder-free nitrile glove, ambidextrous and with rolled rim for extra cuff strength
- Good viral barrier properties, protects against variety of chemicals
- Excellent tactile sensitivity, enhanced wet and dry grip with microtextured fingertips
- Excellent fit, easy to don and remove
- Free of latex, powder and proteins, good alternative to those suffering from type I allergies
- The gloves are suitable for all food contact



Color:  
Black

### Recommendations for use:



Basic Procedures  
and wards



Cleaning



Food Handling

MEDICAL	2017/745	CLASS I	EN 455	AQL 1.5
PPE	2016/425	CATEGORY III	EN 374	AQL 1.5

## MAIN INFORMATION

SIZE	REF	PALM WIDTH	OVERALL LENGTH	THICKNESS		
				PALM	CUFF	FINGER
XS	NIGLXS100BK	75 ± 5 mm	≥240 mm	Min 0.08 mm	Min 0.07 mm	Min 0.10 mm
S	NIGLS100BK	85 ± 5 mm				
M	NIGLM100BK	95 ± 5 mm				
L	NIGLL100BK	110 ± 5 mm				
XL	NIGLXL100BK	120 ± 5 mm				
XXL	NIGLXXL100BK	130 ± 5 mm				

MATERIAL	EXTERIOR	INSIDE TREATMENT	FORCE AT BREAK		SHELF TIME	NATURAL LATEX FREE	ACCELERATOR FREE
			Before ageing	After ageing			
Nitrile	Fingertip textured	Powder Free	≥6N		5 years	Yes	No

## LABELLING

- Each unit displays the following information: product reference, batch number and expiry date.
- Symbols used in compliance with standard EN ISO 15223: Symbols to be used with medical device labels, labelling and information to be supplied.

**CE 2777**



- Notified Body responsible for certification and ongoing conformity: SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland
- Resistance to bacteria and fungi - pass. Resistance to virus - pass.\*The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

## PACKAGING / LOGISTIC INFORMATION

AVAILABLE SIZES	PCS/DISPENSER	PCS/TRANSPORT CARTON	DISPENSER	TRANSPORT CARTON
XS	100 pcs.	1000 pcs.		
S	100 pcs.	1000 pcs.		
M	100 pcs.	1000 pcs.		
L	100 pcs.	1000 pcs.		
XL	100 pcs.	1000 pcs.		
XXL	90 pcs.	900 pcs.		



Sales Unit: Transport Carton

Storage conditions:




**UAB Kasko Group | Kid-Man®**  
 A.Vivulskio str. 12D, 03221, Vilnius, Lithuania  
 www.kid-man.com | info@kid-man.com  
 Tel. +370 5 265 0000

Importator / Distribuitor RO : Kid-Man Aid SRL  
 Strada Theodor Aman Nr 2-2A, Braila Tel: +40 753 660066

## MAIN FEATURES

- Thin powder-free nitrile glove, ambidextrous and with rolled rim for extra cuff strength
- Good viral barrier properties, protects against variety of chemicals
- Excellent tactile sensitivity, enhanced wet and dry grip with microtextured fingertips
- Excellent fit, easy to don and remove
- Free of latex, powder and proteins, good alternative to those suffering from type I allergies
- The gloves are suitable for all food contact



Color:  
Blue

### Recommendations for use:



Basic Procedures and wards



Cleaning



Food Handling

MEDICAL	2017/745	CLASS I	EN 455	AQL 1.5
PPE	2016/425	CATEGORY III	EN 374	AQL 1.5

## MAIN INFORMATION

SIZE	REF	PALM WIDTH	OVERALL LENGTH	THICKNESS		
				PALM	CUFF	FINGER
XS	NIGLXS100BL	75 ± 5 mm	≥240 mm	Min 0.08 mm	Min 0.07 mm	Min 0.10 mm
S	NIGLS100BL	85 ± 5 mm				
M	NIGLM100BL	95 ± 5 mm				
L	NIGLL100BL	110 ± 5 mm				
XL	NIGLXL100BL	120 ± 5 mm				

MATERIAL	EXTERIOR	INSIDE TREATMENT	FORCE AT BREAK		SHELF TIME	NATURAL LATEX FREE	ACCELERATOR FREE
			Before ageing	After ageing			
Nitrile	Fingertip textured	Powder Free	≥6N		5 years	Yes	No

## LABELLING

- Each unit displays the following information: product reference, batch number and expiry date.
- Symbols used in compliance with standard EN ISO 15223: Symbols to be used with medical device labels, labelling and information to be supplied.

**CE 2777**

- Notified Body responsible for certification and ongoing conformity: SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland



**EN ISO 374-5:2016**  
Protects against viruses

- Resistance to bacteria and fungi - pass. Resistance to virus - pass.\*The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

## PACKAGING / LOGISTIC INFORMATION

AVAILABLE SIZES	PCS/DISPENSER	PCS/TRANSPORT CARTON	DISPENSER	TRANSPORT CARTON
XS	100 pcs.	1000 pcs.		
S	100 pcs.	1000 pcs.		
M	100 pcs.	1000 pcs.		
L	100 pcs.	1000 pcs.		
XL	100 pcs.	1000 pcs.		



50 PCS PER BOX



10 PCS PER BOX



1080 BOXES PER PALETTE

Sales Unit: Transport Carton

Storage conditions:




**UAB Kasko Group | Kid-Man®**

A.Vivulskio str. 12D, 03221, Vilnius, Lithuania  
www.kid-man.com | info@kid-man.com  
Tel. +370 5 265 0000

Importator / Distributor RO : Kid-Man Aid SRL  
Strada Theodor Aman Nr 2-2A, Braila Tel: +40 753 660066

## MAIN FEATURES

- Thin powder-free nitrile glove, ambidextrous and with rolled rim for extra cuff strength
- Good viral barrier properties, protects against variety of chemicals
- Excellent tactile sensitivity, enhanced wet and dry grip with microtextured fingertips
- Excellent fit, easy to don and remove
- Free of latex, powder and proteins, good alternative to those suffering from type I allergies
- The gloves are suitable for all food contact



Color:  
Pink

### Recommendations for use:



Basic Procedures and wards



Cleaning



Food Handling

MEDICAL	2017/745	CLASS I	EN 455	AQL 1.5
PPE	2016/425	CATEGORY III	EN 374	AQL 1.5

## MAIN INFORMATION

SIZE	REF	PALM WIDTH	OVERALL LENGTH	THICKNESS		
				PALM	CUFF	FINGER
XS	NIGLXS100RO	75 ± 5 mm	≥240 mm	Min 0.08 mm	Min 0.07 mm	Min 0.10 mm
S	NIGLS100RO	85 ± 5 mm				
M	NIGLM100RO	95 ± 5 mm				
L	NIGLL100RO	110 ± 5 mm				
XL	NIGLXL100RO	120 ± 5 mm				

MATERIAL	EXTERIOR	INSIDE TREATMENT	FORCE AT BREAK		SHELF TIME	NATURAL LATEX FREE	ACCELERATOR FREE
			Before ageing	After ageing			
Nitrile	Fingertip textured	Powder Free	≥6N		5 years	Yes	No

## LABELLING

- Each unit displays the following information: product reference, batch number and expiry date.
- Symbols used in compliance with standard EN ISO 15223: Symbols to be used with medical device labels, labelling and information to be supplied.

**CE 2777**

- Notified Body responsible for certification and ongoing conformity: SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland



**EN ISO 374-5:2016**  
Protects against viruses

- Resistance to bacteria and fungi - pass. Resistance to virus - pass.\*The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

## PACKAGING / LOGISTIC INFORMATION

AVAILABLE SIZES	PCS/DISPENSER	PCS/TRANSPORT CARTON	DISPENSER	TRANSPORT CARTON
XS	100 pcs.	1000 pcs.		
S	100 pcs.	1000 pcs.		
M	100 pcs.	1000 pcs.		
L	100 pcs.	1000 pcs.		
XL	100 pcs.	1000 pcs.		



100 OR 50 PCS PER BOX



10 PCS PER BOX



1080 BOXES PER PALLETTE

Sales Unit: Transport Carton

Storage conditions:



**UAB Kasko Group | Kid-Man®**

A.Vivulskio str. 12D, 03221, Vilnius, Lithuania  
www.kid-man.com | info@kid-man.com  
Tel. +370 5 265 0000

Importator / Distributor RO : Kid-Man Aid SRL  
Strada Theodor Aman Nr 2-2A, Braila Tel: +40 753 660066

## MAIN FEATURES

- Sulphur-free, chemical accelerator-free, latex-free and powder-free
- Skin friendly and reduces the risk of allergies
- Beaded cuff ensures easy donning and prevents rolling down
- Comfortable and reliable hand protection



Color:  
Blue

### Recommendations for use:



Basic Procedures  
and wards



Cleaning



Food Handling

MEDICAL	2017/745	CLASS I	EN 455	AQL 1.5
PPE	2016/425	CATEGORY III	EN 374	AQL 1.5

## MAIN INFORMATION

SIZE	REF	PALM WIDTH	OVERALL LENGTH	THICKNESS		
				PALM	CUFF	FINGER
XS	NIGLXS100AF	≤ 80 mm	≥240 mm	0,07 ± 0,03 mm	0,06 ± 0,03 mm	0,10 ± 0,03 mm
S	NIGLS100AF	80 ±10 mm				
M	NIGLM100AF	95 ±10 mm				
L	NIGLL100AF	110 ±10 mm				
XL	NIGLXL100AF	120 ±10 mm				

MATERIAL	EXTERIOR	INSIDE TREATMENT	FORCE AT BREAK		SHELF TIME	NATURAL LATEX FREE	ACCELERATOR FREE
			Before ageing	After ageing			
Nitrile	Fingertip textured	Chlorinated	≥6N		5 years	Yes	Yes

## LABELLING

- Each unit displays the following information: product reference, batch number and expiry date.
- Symbols used in compliance with standard EN ISO 15223: Symbols to be used with medical device labels, labelling and information to be supplied.

**CE 2777**

- Notified Body responsible for certification and ongoing conformity: SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland



**EN ISO 374-5:2016**  
Protects against viruses

- Resistance to bacteria and fungi - pass. Resistance to virus - pass.\*The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

## PACKAGING / LOGISTIC INFORMATION

AVAILABLE SIZES	PCS/DISPENSER	PCS/TRANSPORT CARTON	DISPENSER	TRANSPORT CARTON
XS	100 pcs.	1000 pcs.		
S	100 pcs.	1000 pcs.		
M	100 pcs.	1000 pcs.		
L	100 pcs.	1000 pcs.		
XL	100 pcs.	1000 pcs.		



100 OR 50 PCS PER BOX



10 PCS PER BOX



800 BOXES PER PALLETTE

Sales Unit: Transport Carton

Storage conditions:



**UAB Kasko Group | Kid-Man®**

A.Vivulskio str. 12D, 03221, Vilnius, Lithuania  
www.kid-man.com | info@kid-man.com  
Tel. +370 5 265 0000

Importator / Distributor RO : Kid-Man Aid SRL  
Strada Theodor Aman Nr 2-2A, Braila Tel:Tel: +40 753 660066

## MAIN FEATURES

- Extended cuff to offer better protection by enclosing the sleeve
- Excellent softness and wearing fitness
- Finger textured for a reliable, consistent grip
- Beaded cuff ensures easy donning and prevents roll down
- Latex-free, powder-free and odor-free
- Compliant with food handling requirements



Color:  
Blue

### Recommendations for use:



Basic Procedures and wards



Cleaning



Food Handling

MEDICAL	2017/745	CLASS I	EN 455	AQL 1.5
PPE	2016/425	CATEGORY III	EN 374	AQL 1.5

## MAIN INFORMATION

SIZE	REF	PALM WIDTH	OVERALL LENGTH	THICKNESS		
				PALM	CUFF	FINGER
XS	NIGLXS100290	70 ±10 mm	≥ 290 mm	0,09 ± 0,03 mm	0,08 ± 0,03 mm	0,14 ± 0,03 mm
S	NIGLS100290	80 ±10 mm				
M	NIGLM100290	95 ±10 mm				
L	NIGLL100290	110 ±10 mm				
XL	NIGLXL100290	120 ±10 mm				

MATERIAL	EXTERIOR	INSIDE TREATMENT	FORCE AT BREAK		SHELF TIME	NATURAL LATEX FREE	ACCELERATOR FREE
			Before ageing	After ageing			
Nitrile	Fingertip textured	Chlorinated	≥6N		5 years	Yes	No

## LABELLING

- Each unit displays the following information: product reference, batch number and expiry date.
- Symbols used in compliance with standard EN ISO 15223: Symbols to be used with medical device labels, labelling and information to be supplied.

**CE 2777**

- Notified Body responsible for certification and ongoing conformity: SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland



**EN ISO 374-5:2016**  
Protects against viruses

- Resistance to bacteria and fungi - pass. Resistance to virus - pass.\*The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

## PACKAGING / LOGISTIC INFORMATION

AVAILABLE SIZES	PCS/DISPENSER	PCS/TRANSPORT CARTON	DISPENSER	TRANSPORT CARTON
XS	100 pcs.	1000 pcs.		
S	100 pcs.	1000 pcs.		
M	100 pcs.	1000 pcs.		
L	100 pcs.	1000 pcs.		
XL	100 pcs.	1000 pcs.		



100 OR 50 PCS PER BOX



10 PCS PER BOX



800 BOXES PER PALETTE

**Sales Unit:** Transport Carton

**Storage conditions:**



**UAB Kasko Group | Kid-Man®**

A.Vivulskio str. 12D, 03221, Vilnius, Lithuania  
www.kid-man.com | info@kid-man.com  
Tel. +370 5 265 0000

Importator / Distributor RO : Kid-Man Aid SRL  
Strada Theodor Aman Nr 2-2A, Braila Tel: +40 753 660066

# Biodegradable Nitrile Gloves

## MAIN FEATURES

- Biodegradable gloves
- Biodegradation within 270 days
- Contains innovative biodegradable technology
- Powder free



Color:  
Blue

Type:  
Biodegradable

Recommendations for use:



Basic Procedures  
and wards



Cleaning



Food Handling

MEDICAL	2017/745	CLASS I	EN 455	AQL 1.5
PPE	2016/425	CATEGORY III	EN 374	AQL 1.5

## MAIN INFORMATION

SIZE	REF	PALM WIDTH	OVERALL LENGTH	THICKNESS		
				PALM	CUFF	FINGER
XS	NIGLXS100BDG	76 ± 4 mm	≥240 mm	0,07 ± 0,02 mm	0,06 ± 0,03 mm	0,10 ± 0,03 mm
S	NIGLS100BDG	86 ± 4 mm				
M	NIGLM100BDG	98 ± 4 mm				
L	NIGLL100BDG	107 ± 4 mm				
XL	NIGLXL100BDG	115 ± 4 mm				

MATERIAL	EXTERIOR	INSIDE TREATMENT	FORCE AT BREAK		SHELF TIME	NATURAL LATEX FREE	ACCELERATOR FREE
			Before ageing	After ageing			
Nitrile Latex	Fingertip textured	Chlorinated	≥6N		5 years	Yes	No

# Biodegradable Nitrile Gloves

## LABELLING

- Each unit displays the following information: product reference, batch number and expiry date.
- Symbols used in compliance with standard EN ISO 15223: Symbols to be used with medical device labels, labelling and information to be supplied.

**CE 2777**


- Notified Body responsible for certification and ongoing conformity: SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland



- Resistance to bacteria and fungi - pass. Resistance to virus - pass.\*The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

## PACKAGING / LOGISTIC INFORMATION

AVAILABLE SIZES	PCS/DISPENSER	PCS/TRANSPORT CARTON	DISPENSER	TRANSPORT CARTON
XS	100 pcs.	1000 pcs.		
S	100 pcs.	1000 pcs.		
M	100 pcs.	1000 pcs.		
L	100 pcs.	1000 pcs.		
XL	100 pcs.	1000 pcs.		



Sales Unit: Transport Carton

Storage conditions:




**UAB Kasko Group | Kid-Man®**  
 A.Vivulskio str. 12D, 03221, Vilnius, Lithuania  
 www.kid-man.com | info@kid-man.com  
 Tel. +370 5 265 0000

Importator / Distribuitor RO : Kid-Man Aid SRL  
 Strada Theodor Aman Nr 2-2A, Braila Tel:Tel: +40 753 660066

KID-MAN GRIP gloves are designed for professionals working in demanding environments. Made from flexible and durable nitrile, these gloves are twice as thick as standard thin nitrile gloves, offering enhanced protection. The unique diamond 3D texture ensures a secure grip, even when handling wet or oily items.


## MAIN FEATURES

- **Strong grip** – the diamond texture provides stability in greasy and wet environments.
- **Comfortable wear** – inner texture reduces hand sweating.
- **High protection** – shields against select chemicals, bacteria, viruses, and fungi.
- **Durability** – less wear and tear than standard nitrile gloves.
- **Food-safe** – approved for food processing applications.



Color:  
Black

### Recommendations for use:



MEDICAL	2017/745	CLASS I	EN 455	AQL 1.5
PPE	2016/425	CATEGORY III	EN 374	AQL 1.5

## MAIN INFORMATION

SIZE	REF	PALM WIDTH	OVERALL LENGTH	THICKNESS		
				PALM	CUFF	FINGER
M	NIGLM50GRIP	95 ±10 mm	≥240 mm	0,17 ± 0,03 mm	0,09 ± 0,03 mm	0,21 ± 0,03 mm
L	NIGLL50GRIP	110 ±10 mm				
XL	NIGLXL50GRIP	120 ±10 mm				
XXL	NIGLXXL50GRIP	130 ±10 mm				

MATERIAL	EXTERIOR	INSIDE TREATMENT	FORCE AT BREAK		SHELF TIME	NATURAL LATEX FREE	ACCELERATOR FREE
			Before ageing	After ageing			
Nitrile	Diamond textured	Chlorinated	≥6N		5 years	Yes	No

## LABELLING

- Each unit displays the following information: product reference, batch number and expiry date.
- Symbols used in compliance with standard EN ISO 15223: Symbols to be used with medical device labels, labelling and information to be supplied.

**CE 2777**

- Notified Body responsible for certification and ongoing conformity: SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland



- Resistance to bacteria and fungi - pass. Resistance to virus - pass.\*The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

## PACKAGING / LOGISTIC INFORMATION

AVAILABLE SIZES	PCS/DISPENSER	PCS/TRANSPORT CARTON	DISPENSER	TRANSPORT CARTON
M	50 pcs.	500 pcs.		
L	50 pcs.	500 pcs.		
XL	50 pcs.	500 pcs.		
XXL	50 pcs.	500 pcs.		



Sales Unit: Transport Carton

Storage conditions:



**UAB Kasko Group | Kid-Man®**  
 A.Vivulskio str. 12D, 03221, Vilnius, Lithuania  
 www.kid-man.com | info@kid-man.com  
 Tel. +370 5 265 0000

Importator / Distribuitor RO : Kid-Man Aid SRL  
 Strada Theodor Aman Nr 2-2A, Braila Tel:Tel: +40 753 660066

## MAIN FEATURES

- Latex gloves are made of natural rubber (NR) latex, which offers good resistance against many acids and alkalis
- The best fit and feel, latex gloves are very comfortable to wear
- A disadvantage of natural latex is that its natural proteins can cause or trigger latex allergies
- Protein content does not exceed 50 µg/g



Color:  
White

### Recommendations for use:



Basic Procedures  
and wards



Cleaning



Food Handling

MEDICAL	2017/745	CLASS I	EN 455	AQL 1.5
PPE	2016/425	CATEGORY III	EN 374	AQL 1.5

## MAIN INFORMATION

SIZE	REF	PALM WIDTH	OVERALL LENGTH	THICKNESS		
				PALM	CUFF	FINGER
XS	LTGLXS100	74 ± 4 mm	≥240 mm	0,10 ± 0,02 mm	0,08 ± 0,03 mm	0,12 ± 0,03 mm
S	LTGLS100	84 ± 4 mm				
M	LTGLM100	94 ± 4 mm				
L	LTGLL100	105 ± 4 mm				
XL	LTGLXL100	114 ± 4 mm				

MATERIAL	EXTERIOR	INSIDE TREATMENT	FORCE AT BREAK		SHELF TIME	NATURAL LATEX FREE	ACCELERATOR FREE
			Before ageing	After ageing			
Natural Rubber Latex	Fingertip textured	Chlorinated	≥6N		5 years	No	No

## LABELLING

- Each unit displays the following information: product reference, batch number and expiry date.
- Symbols used in compliance with standard EN ISO 15223: Symbols to be used with medical device labels, labelling and information to be supplied.

**CE 2777**

- Notified Body responsible for certification and ongoing conformity: SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland



**EN ISO 374-5:2016**  
Protects against viruses

- Resistance to bacteria and fungi - pass. Resistance to virus - pass.\*The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

## PACKAGING / LOGISTIC INFORMATION

AVAILABLE SIZES	PCS/DISPENSER	PCS/TRANSPORT CARTON	DISPENSER	TRANSPORT CARTON
XS	100 pcs.	1000 pcs.		
S	100 pcs.	1000 pcs.		
M	100 pcs.	1000 pcs.		
L	100 pcs.	1000 pcs.		
XL	100 pcs.	1000 pcs.		



100 OR 50 PCS PER BOX



10 PCS PER BOX



800 BOXES PER PALLETTE

**Sales Unit:** Transport Carton

**Storage conditions:**



**UAB Kasko Group | Kid-Man®**  
A.Vivulskio str. 12D, 03221, Vilnius, Lithuania  
www.kid-man.com | info@kid-man.com  
Tel. +370 5 265 0000

Importator / Distributor RO : Kid-Man Aid SRL  
Strada Theodor Aman Nr 2-2A, Braila Tel: +40 753 660066

## MAIN FEATURES

- Thin powder-free vinyl glove, ambidextrous and with rolled rim for extra cuff strength
- Good viral barrier properties, protects against variety of chemicals
- Easy to don and remove
- Free of latex, powder and proteins, good alternative to those suffering from type I allergies
- The gloves are suitable for all food contact



Color:  
Clear

### Recommendations for use:



Basic Procedures and wards



Cleaning



Food Handling

MEDICAL	2017/745	CLASS I	EN 455	AQL 1.5
PPE	2016/425	CATEGORY III	EN 374	AQL 1.5

## MAIN INFORMATION

SIZE	REF	PALM WIDTH	OVERALL LENGTH	THICKNESS		
				PALM	CUFF	FINGER
S	VIGLS100	85 ± 5 mm	≥240 mm	Min 0.08 mm	Min 0.06 mm	Min 0.10 mm
M	VIGLM100	95 ± 5 mm				
L	VIGLL100	105 ± 5 mm				
XL	VIGLXL100	115 ± 5 mm				

MATERIAL	EXTERIOR	INSIDE TREATMENT	FORCE AT BREAK		SHELF TIME	NATURAL LATEX FREE	ACCELERATOR FREE
			Before ageing	After ageing			
Polyvinyl chloride (PVC)	Smooth	Powder Free	Min. 3.6N		5 years	Yes	No

## LABELLING

- Each unit displays the following information: product reference, batch number and expiry date.
- Symbols used in compliance with standard EN ISO 15223: Symbols to be used with medical device labels, labelling and information to be supplied.

**CE 2777**

- Notified Body responsible for certification and ongoing conformity: SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland



- Resistance to bacteria and fungi - pass. Resistance to virus - pass.\*The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

## PACKAGING / LOGISTIC INFORMATION

AVAILABLE SIZES	PCS/DISPENSER	PCS/TRANSPORT CARTON	DISPENSER	TRANSPORT CARTON
S	100 pcs.	1000 pcs.		
M	100 pcs.	1000 pcs.		
L	100 pcs.	1000 pcs.		
XL	100 pcs.	1000 pcs.		



Sales Unit: Transport Carton

Storage conditions:



**UAB Kasko Group | Kid-Man®**  
 A.Vivulskio str. 12D, 03221, Vilnius, Lithuania  
 www.kid-man.com | info@kid-man.com  
 Tel. +370 5 265 0000

Importator / Distribuitor RO : Kid-Man Aid SRL  
 Strada Theodor Aman Nr 2-2A, Braila Tel:Tel: +40 753 660066


## Regulatory & requirements

What do the pictograms and symbols on the glove boxes mean?

 Manufacturer	 Date of manufacture	 Temperature limitation	 Consult instructions for use
 Lot / batch number	 Expiry date	 For single use only	 Recycling packaging
 Product number	 Keep dry	 Non-sterile	 This side up
 Medical device	 Keep away from sunlight	 Suitable for food contact	 Handle with care
 Serial number			

## Let's get in touch

 UAB Kasko Group | Kid-Man®  
A.Vivulskio str. 12D  
03221, Vilnius, Lithuania  
Tel. +370 5 265 0000  
order@kid-man.com

 Kid-Man AID S.R.L  
Str. Theodor Aman nr.2-2A  
810031 Braila, Romania  
Telefon +40 753 660066  
ro.order@kid-man.com



[www.kid-man.com](http://www.kid-man.com)