

Kid-Man® **Examination gloves**

TECHNICAL DATA SHEET





Key features of NITRILE GLOVES

- Thin powder-free nitrile glove, ambidextrous and with rolled rim for extra cuff strength
- Good viral barrier properties, protects against a variety of chemicals
- Excellent tactile sensitivity, enhanced wet and dry grip with microtextured fingertips
- Extremely tear-resistant
- Excellent fit, easy to don and remove
- Free of latex, powder and proteins, good alternative to those suffering from type I allergies
- The gloves are suitable for all food contact

NIGL100RO (XS-XL)



NIGL100BL (XS-XL)



NIGL100 (XS-XL)

NIGL100BK (XS-XXL)

MEDICAL	2017/745	CLASS I	EN 455	AQL 1.5
PPE	2016/425	CATEGORY III	EN 374	AQL 1.5



NIGL100AF (XS-XL)

Accelerator-free For allergic skin

- Sulphur-free, chemical accelerator-free, latex-free and powder-free
- Skin friendly and reduces the risk of allergies
- Beaded cuff ensures easy donning and prevents rolling down
- Comfortable and reliable hand protection

Application areas

- Chemical industry
- Workshop activities
- Painting, varnishing and polishing work
- Electronics and computer industry
- Laboratory and research
- Cleaning and maintenance
- Product protection and printing industry



NIGL100290 (XS-XL)

290 mm length More safety & comfort

- Longer and stronger
- Sulphur-free, latex-free and powder-free
- Skin friendly and reduces the risk of allergies
- Beaded cuff ensures easy donning and prevents rolling down
- Comfortable and reliable hand protection

Application areas

- Chemical industry
- Workshop activities
- Painting, varnishing and polishing work
- Electronics and computer industry
- Laboratory and research
- Cleaning and maintenance
- Product protection and printing industry



NIGL50GRIP (M-XXL)

Diamond texture For optimal grip

- Extremely tear-resistant
- Silicone-free manufacturing process
- Reliable hand protection

Application areas

- Heavy-duty work, vehicle repair
- Chemical industry
- Workshop activities
- Painting, varnishing and polishing work
- Electronics and computer industry
- Laboratory and research
- Cleaning and maintenance
- Product protection and printing industry



NIGL100BDG (XS-XL)

Biodegradable Eco-friendly choice

- Biodegradable gloves
- Biodegradation within 270 days
- Contains innovative biodegradable technology
- Comfortable and reliable hand protection

Application areas

- Green and sustainable companies
- Chemical industry
- Workshop activities
- Painting, varnishing and polishing work
- Electronics and computer industry
- Laboratory and research
- Cleaning and maintenance
- Product protection and printing industry

Key features of VINYL & LATEX gloves



LTGL100 (XS-XL)

Latex

- Made of natural rubber (NR) latex, which offers good resistance against many acids and alkalis
- The best fit and feel, comfortable to wear
- A disadvantage of natural latex is that its natural proteins can cause or trigger latex allergies.
- Protein content does not exceed 50 µg/g.

Application areas

- Healthcare, nursing and care homes
- Dental clinics
- Beauty and tattoo studios
- Food industry and catering
- Laboratory and research
- Cleaning and maintenance
- Household and janitorial work



VIGL100 (S-XL)

Vinyl

- Thin powder-free glove, ambidextrous and with rolled rim for extra cuff strength
- Good viral barrier properties, protects against a variety of chemicals
- Easy to don and remove
- Free of latex, powder and proteins, good alternative to those suffering from type I allergies
- The gloves are suitable for all food contact.

Application areas

- Food industry and catering (popular among restaurant, takeaway and sushi bar kitchens)
- Beauty and hair salons
- Cleaning and maintenance
- Retail and cashiers
- Light workshop activities
- Household and janitorial work
- Product protection and packing

Technical details about Kid-Man® gloves assortment

Size range

What sizes of gloves can you get?



Colour range

What colours of gloves can you get?



Which glove material should I choose?

Individuals may perceive the material properties of nitrile, vinyl and latex differently. The following comparison provides a general overview and summarises the main material characteristics to help you choose the right type of glove.

Does it match what you are looking for?  The best choice  Recommended  Not recommended

Material	NITRILE								VINYL	LATEX
REF / Product name	NIGL100 NIGL100BL NIGL100RO	NIGL100BK	NIGL100290 NIGL100LIGHT	NIGL100AF	NIGL100BDG	NIGL100OAT	NIGL100EKO	NIGL50GRIP	VIGL100	LTGL100
Comfort	✓	✓	✓	✓	✓	✓	✓	✓	✗	
Elasticity	✓	✓	✓	✓	✓	✓	✓	✓	✗	✓
Grip / Tactility	✓	✓	✓	✓	✓	✓	✓		✗	✓
Tear resistance	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Elongation	600 %	600 %	600 %	600 %	600 %	600 %	600 %	600 %	300 %	800 %
Puncture resistance	✓	✓	✓	✓	✓	✓	✓		✓	✓
Chemical resistance	✓	✓	✓	✓	✓	✓	✓	✓	✗	✓
Latex-free	✓	✓	✓	✓	✓	✓	✓	✓	✓	✗
Accelerator-free	✗	✗	✗		✗	✗	✗	✗	✓	✗
Plasticisers	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Heavy duty / extra thick	✗	✗	✗	✗	✗	✗	✗		✓	✗
Ecological	✗	✗	✗	✗		✗	✗	✗	✗	✗
Size range	XS-XL	XS-XXL	XS-XL	XS-XL	XS-XL	XS-XL	XS-XXL	M-XXL	S-XL	XS-XL
Colour	Violet Blue Pink	Black	Blue D. Blue	Blue	D. Blue	Blue	Blue	Black	Clear	White

What are the logistics of transporting the gloves?

Sales unit: transport carton.



REF / MODEL NAME	GLOVES PER PACKAGE (PCS.)	PACKAGES PER TRANSPORT CARTON (PCS.)	PACKAGES ON THE PALLET (PCS.)
NIGL100, NIGL100BK, NIGL100BL, NIGL100RO, NIGL100LIGHT	100 (XXL size 90 pcs.)	10	1080
NIGL100AF, VIGL, LTGL	100	10	800
NIGL100290	100	10	630
NIGL100BDG	100	10	840
NIGL50GRIP	50	10	800

MAIN FEATURES

- Thin powder-free nitrile glove, ambidextrous and with rolled rim for extra cuff strength
- Good viral barrier properties, protects against variety of chemicals
- Excellent tactile sensitivity, enhanced wet and dry grip with microtextured fingertips
- Excellent fit, easy to don and remove
- Free of latex, powder and proteins, good alternative to those suffering from type I allergies
- The gloves are suitable for all food contact



Color:
Violet

Recommendations for use:



Basic Procedures
and wards



Cleaning



Food Handling

MEDICAL	2017/745	CLASS I	EN 455	AQL 1.5
PPE	2016/425	CATEGORY III	EN 374	AQL 1.5

MAIN INFORMATION

SIZE	REF	PALM WIDTH	OVERALL LENGHT	THICKNESS			
				PALM	CUFF	FINGER	
XS	NIGLXS100	75 ± 5 mm	≥240 mm	Min 0.08 mm	Min 0.07 mm	Min 0.10 mm	
S	NIGLS100	85 ± 5 mm					
M	NIGLM100	95 ± 5 mm					
L	NIGLL100	110 ± 5 mm					
XL	NIGLXL100	120 ± 5 mm					
MATERIAL	EXTERIOR	INSIDE TREATMENT	FORCE AT BREAK		SHELF TIME	NATURAL LATEX FREE	ACCELERATOR FREE
			Before ageing	After ageing			

LABELLING

- Each unit displays the following information: product reference, batch number and expiry date.
- Symbols used in compliance with standard EN ISO 15223: Symbols to be used with medical device labels, labelling and information to be supplied.

CE 2777

- Notified Body responsible for certification and ongoing conformity: SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland



EN ISO 374-5:2016
Protects against viruses

- Resistance to bacteria and fungi - pass. Resistance to virus - pass.*The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

PACKAGING / LOGISTIC INFORMATION

AVAILABLE SIZES	PCS/DISPENSER	PCS/TRANSPORT CARTON	DISPENSER	TRANSPORT CARTON
XS	100 pcs.	1000 pcs.	 4 779031 632806	 4 779031 632738
S	100 pcs.	1000 pcs.	 4 779031 631984	 4 779031 632721
M	100 pcs.	1000 pcs.	 4 779031 631960	 4 779031 632752
L	100 pcs.	1000 pcs.	 4 779031 631977	 4 779031 632776
XL	100 pcs.	1000 pcs.	 4 779031 631953	 4 779031 632745



100 PCS PER BOX



10 PCS PER BOX



1080 BOXES PER PALLETTE

Sales Unit: Transport Carton

Storage conditions:




UAB Kasko Group | Kid-Man®

A.Vivulskio str. 12D, 03221, Vilnius, Lithuania
www.kid-man.com | info@kid-man.com
Tel. +370 5 265 0000

Importator / Distribuitor RO : Kid-Man Aid SRL
Strada Theodor Aman Nr 2-2A, Braila Tel: +40 753 660066

MAIN FEATURES

- Thin powder-free nitrile glove, ambidextrous and with rolled rim for extra cuff strength
- Good viral barrier properties, protects against variety of chemicals
- Excellent tactile sensitivity, enhanced wet and dry grip with microtextured fingertips
- Excellent fit, easy to don and remove
- Free of latex, powder and proteins, good alternative to those suffering from type I allergies
- The gloves are suitable for all food contact



Color:
Black

Recommendations for use:



Basic Procedures
and wards



Cleaning



Food Handling

MEDICAL	2017/745	CLASS I	EN 455	AQL 1.5
PPE	2016/425	CATEGORY III	EN 374	AQL 1.5

MAIN INFORMATION

SIZE	REF	PALM WIDTH	OVERALL LENGHT	THICKNESS			
				PALM	CUFF	FINGER	
XS	NIGLXS100BK	75 ± 5 mm	≥240 mm	Min 0.08 mm	Min 0.07 mm	Min 0.10 mm	
S	NIGLS100BK	85 ± 5 mm					
M	NIGLM100BK	95 ± 5 mm					
L	NIGLL100BK	110 ± 5 mm					
XL	NIGLXL100BK	120 ± 5 mm					
XXL	NIGLXXL100BK	130 ± 5 mm					
MATERIAL	EXTERIOR	INSIDE TREATMENT	FORCE AT BREAK		SHELF TIME	NATURAL LATEX FREE	ACCELERATOR FREE
			Before ageing	After ageing			
Nitrile	Fingertip textured	Powder Free	≥6N	≥6N	5 years	Yes	No

LABELLING

- Each unit displays the following information: product reference, batch number and expiry date.
- Symbols used in compliance with standard EN ISO 15223: Symbols to be used with medical device labels, labelling and information to be supplied.

CE 2777



- Notified Body responsible for certification and ongoing conformity: SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland
- Resistance to bacteria and fungi - pass. Resistance to virus - pass.*The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

PACKAGING / LOGISTIC INFORMATION

AVAILABLE SIZES	PCS/DISPENSER	PCS/TRANSPORT CARTON	DISPENSER	TRANSPORT CARTON
XS	100 pcs.	1000 pcs.	 4 779031 632820	 4 779031 632875
S	100 pcs.	1000 pcs.	 4 779031 632813	 4 779031 632868
M	100 pcs.	1000 pcs.	 4 779031 632837	 4 779031 632882
L	100 pcs.	1000 pcs.	 4 779031 632844	 4 779031 632899
XL	100 pcs.	1000 pcs.	 4 779031 632851	 4 779031 632905
XXL	90 pcs.	900 pcs.	 4 779031 633254	 4 779031 633261



100 PCS PER BOX



10 PCS PER BOX



1080 BOXES PER PALLETTE

Sales Unit: Transport Carton

Storage conditions:



UAB Kasko Group | Kid-Man®

A.Vivulskio str. 12D, 03221, Vilnius, Lithuania
www.kid-man.com | info@kid-man.com
Tel. +370 5 265 0000

Importator / Distribuitor RO : Kid-Man Aid SRL
Strada Theodor Aman Nr 2-2A, Braila Tel: +40 753 660066

MAIN FEATURES

- Thin powder-free nitrile glove, ambidextrous and with rolled rim for extra cuff strength
- Good viral barrier properties, protects against variety of chemicals
- Excellent tactile sensitivity, enhanced wet and dry grip with microtextured fingertips
- Excellent fit, easy to don and remove
- Free of latex, powder and proteins, good alternative to those suffering from type I allergies
- The gloves are suitable for all food contact



Color:
Blue

Recommendations for use:



Basic Procedures
and wards



Cleaning



Food Handling

MEDICAL	2017/745	CLASS I	EN 455	AQL 1.5
PPE	2016/425	CATEGORY III	EN 374	AQL 1.5

MAIN INFORMATION

SIZE	REF	PALM WIDTH	OVERALL LENGHT	THICKNESS			
				PALM	CUFF	FINGER	
XS	NIGLXS100BL	75 ± 5 mm	≥240 mm	Min 0.08 mm	Min 0.07 mm	Min 0.10 mm	
S	NIGLS100BL	85 ± 5 mm					
M	NIGLM100BL	95 ± 5 mm					
L	NIGLL100BL	110 ± 5 mm					
XL	NIGLXL100BL	120 ± 5 mm					
MATERIAL	EXTERIOR	INSIDE TREATMENT	FORCE AT BREAK		SHELF TIME	NATURAL LATEX FREE	ACCELERATOR FREE
			Before ageing	After ageing			
Nitrile	Fingertip textured	Powder Free	≥6N		5 years	Yes	No

LABELLING

- Each unit displays the following information: product reference, batch number and expiry date.
- Symbols used in compliance with standard EN ISO 15223: Symbols to be used with medical device labels, labelling and information to be supplied.

CE 2777

- Notified Body responsible for certification and ongoing conformity: SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland



EN ISO 374-5:2016
Protects against viruses

- Resistance to bacteria and fungi - pass. Resistance to virus - pass.*The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

PACKAGING / LOGISTIC INFORMATION

AVAILABLE SIZES	PCS/DISPENSER	PCS/TRANSPORT CARTON	DISPENSER	TRANSPORT CARTON
XS	100 pcs.	1000 pcs.	 4 779031 632936	 4 779031 633032
S	100 pcs.	1000 pcs.	 4 779031 632943	 4 779031 633049
M	100 pcs.	1000 pcs.	 4 779031 632950	 4 779031 633056
L	100 pcs.	1000 pcs.	 4 779031 632967	 4 779031 633063
XL	100 pcs.	1000 pcs.	 4 779031 632974	 4 779031 633070



100 PCS PER BOX



10 PCS PER BOX



1080 BOXES PER PALLETTE

Sales Unit: Transport Carton

Storage conditions:



UAB Kasko Group | Kid-Man®

A.Vivulskio str. 12D, 03221, Vilnius, Lithuania
www.kid-man.com | info@kid-man.com
Tel. +370 5 265 0000

Importator / Distribuitor RO : Kid-Man Aid SRL
Strada Theodor Aman Nr 2-2A, Braila Tel: +40 753 660066

MAIN FEATURES

- Thin powder-free nitrile glove, ambidextrous and with rolled rim for extra cuff strength
- Good viral barrier properties, protects against variety of chemicals
- Excellent tactile sensitivity, enhanced wet and dry grip with microtextured fingertips
- Excellent fit, easy to don and remove
- Free of latex, powder and proteins, good alternative to those suffering from type I allergies
- The gloves are suitable for all food contact



Color:
Pink

Recommendations for use:



Basic Procedures
and wards



Cleaning



Food Handling

MEDICAL	2017/745	CLASS I	EN 455	AQL 1.5
PPE	2016/425	CATEGORY III	EN 374	AQL 1.5

MAIN INFORMATION

SIZE	REF	PALM WIDTH	OVERALL LENGHT	THICKNESS			
				PALM	CUFF	FINGER	
XS	NIGLXS100RO	75 ± 5 mm	≥240 mm	Min 0.08 mm	Min 0.07 mm	Min 0.10 mm	
S	NIGLS100RO	85 ± 5 mm					
M	NIGLM100RO	95 ± 5 mm					
L	NIGLL100RO	110 ± 5 mm					
XL	NIGLXL100RO	120 ± 5 mm					
MATERIAL	EXTERIOR	INSIDE TREATMENT	FORCE AT BREAK		SHELF TIME	NATURAL LATEX FREE	ACCELERATOR FREE
			Before ageing	After ageing			
Nitrile	Fingertip textured	Powder Free	≥6N		5 years	Yes	No

LABELLING

- Each unit displays the following information: product reference, batch number and expiry date.
- Symbols used in compliance with standard EN ISO 15223: Symbols to be used with medical device labels, labelling and information to be supplied.

CE 2777

- Notified Body responsible for certification and ongoing conformity: SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland



EN ISO 374-5:2016
Protects against viruses

- Resistance to bacteria and fungi - pass. Resistance to virus - pass.*The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

PACKAGING / LOGISTIC INFORMATION

AVAILABLE SIZES	PCS/DISPENSER	PCS/TRANSPORT CARTON	DISPENSER	TRANSPORT CARTON
XS	100 pcs.	1000 pcs.	 4 779031 634510	 4 779031 634572
S	100 pcs.	1000 pcs.	 4 779031 634527	 4 779031 634589
M	100 pcs.	1000 pcs.	 4 779031 634534	 4 779031 634596
L	100 pcs.	1000 pcs.	 4 779031 634541	 4 779031 634602
XL	100 pcs.	1000 pcs.	 4 779031 634558	 4 779031 634619



100 PCS PER BOX



10 PCS PER BOX



1080 BOXES PER PALLETTE

Sales Unit: Transport Carton

Storage conditions:



UAB Kasko Group | Kid-Man®

A.Vivulskio str. 12D, 03221, Vilnius, Lithuania
www.kid-man.com | info@kid-man.com
Tel. +370 5 265 0000

Importator / Distribuitor RO : Kid-Man Aid SRL
Strada Theodor Aman Nr 2-2A, Braila Tel: +40 753 660066

MAIN FEATURES

- Sulphur-free, chemical accelerator-free, latex-free and powder-free
- Skin friendly and reduces the risk of allergies
- Beaded cuff ensures easy donning and prevents rolling down
- Comfortable and reliable hand protection



Color:
Blue

Recommendations for use:



Basic Procedures
and wards



Cleaning



Food Handling

MEDICAL	2017/745	CLASS I	EN 455	AQL 1.5
PPE	2016/425	CATEGORY III	EN 374	AQL 1.5

MAIN INFORMATION

SIZE	REF	PALM WIDTH	OVERALL LENGHT	THICKNESS		
				PALM	CUFF	FINGER
XS	NIGLXS100AF	≤ 80 mm	≥240 mm	0.07 ± 0,03 mm	0.06 ± 0,03 mm	0.10 ± 0.03 mm
S	NIGLS100AF	80 ±10 mm				
M	NIGLM100AF	95 ±10 mm				
L	NIGLL100AF	110 ±10 mm				
XL	NIGLXL100AF	120 ±10 mm				

MATERIAL	EXTERIOR	INSIDE TREATMENT	FORCE AT BREAK		SHELF TIME	NATURAL LATEX FREE	ACCELERATOR FREE
			Before ageing	After ageing			
Nitrile	Fingertip textured	Chlorinated	≥6N		5 years	Yes	Yes

LABELLING

- Each unit displays the following information: product reference, batch number and expiry date.
- Symbols used in compliance with standard EN ISO 15223: Symbols to be used with medical device labels, labelling and information to be supplied.

CE 2777

- Notified Body responsible for certification and ongoing conformity: SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland



EN ISO 374-5:2016
Protects against viruses

- Resistance to bacteria and fungi - pass. Resistance to virus - pass.*The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

PACKAGING / LOGISTIC INFORMATION

AVAILABLE SIZES	PCS/DISPENSER	PCS/TRANSPORT CARTON	DISPENSER	TRANSPORT CARTON
XS	100 pcs.	1000 pcs.	 4 779031 633681	 4 779031 633674
S	100 pcs.	1000 pcs.	 4 779031 633704	 4 779031 633698
M	100 pcs.	1000 pcs.	 4 779031 633711	 4 779031 633728
L	100 pcs.	1000 pcs.	 4 779031 633742	 4 779031 633735
XL	100 pcs.	1000 pcs.	 4 779031 633766	 4 779031 633759



100 PCS PER BOX



10 PCS PER BOX



800 BOXES PER PALLETTE

Sales Unit: Transport Carton

Storage conditions:



UAB Kasko Group | Kid-Man®

A.Vivulskio str. 12D, 03221, Vilnius, Lithuania
www.kid-man.com | info@kid-man.com
Tel. +370 5 265 0000

Importator / Distribuitor RO : Kid-Man Aid SRL
Strada Theodor Aman Nr 2-2A, Braila Tel: +40 753 660066

MAIN FEATURES

- Extended cuff to offer better protection by enclosing the sleeve
- Excellent softness and wearing fitness
- Finger textured for a reliable, consistent grip
- Beaded cuff ensures easy donning and prevents roll down
- Latex-free, powder-free and odor-free
- Compliant with food handling requirements



Color:
Blue

Recommendations for use:



Basic Procedures
and wards



Cleaning



Food Handling

MEDICAL	2017/745	CLASS I	EN 455	AQL 1.5
PPE	2016/425	CATEGORY III	EN 374	AQL 1.5

MAIN INFORMATION

SIZE	REF	PALM WIDTH	OVERALL LENGHT	THICKNESS		
				PALM	CUFF	FINGER
XS	NIGLXS100290	70 ±10 mm	≥ 290 mm	0.09 ± 0,03 mm	0.08 ± 0,03 mm	0.14 ± 0.03 mm
S	NIGLS100290	80 ±10 mm				
M	NIGLM100290	95 ±10 mm				
L	NIGLL100290	110 ±10 mm				
XL	NIGLXL100290	120 ±10 mm				

MATERIAL	EXTERIOR	INSIDE TREATMENT	FORCE AT BREAK		SHELF TIME	NATURAL LATEX FREE	ACCELERATOR FREE
			Before ageing	After ageing			
Nitrile	Fingertip textured	Chlorinated	≥6N		5 years	Yes	No

LABELLING

- Each unit displays the following information: product reference, batch number and expiry date.
- Symbols used in compliance with standard EN ISO 15223: Symbols to be used with medical device labels, labelling and information to be supplied.

CE 2777

- Notified Body responsible for certification and ongoing conformity: SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland



EN ISO 374-5:2016
Protects against viruses

- Resistance to bacteria and fungi - pass. Resistance to virus - pass.*The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

PACKAGING / LOGISTIC INFORMATION

AVAILABLE SIZES	PCS/DISPENSER	PCS/TRANSPORT CARTON	DISPENSER	TRANSPORT CARTON
XS	100 pcs.	1000 pcs.	 4 779031 633889	 4 779031 633872
S	100 pcs.	1000 pcs.	 4 779031 633902	 4 779031 633896
M	100 pcs.	1000 pcs.	 4 779031 633926	 4 779031 633919
L	100 pcs.	1000 pcs.	 4 779031 633933	 4 779031 633940
XL	100 pcs.	1000 pcs.	 4 779031 633957	 4 779031 633964



100 PCS PER BOX



10 PCS PER BOX



630 BOXES PER PALLETTE

Sales Unit: Transport Carton

Storage conditions:



UAB Kasko Group | Kid-Man®

A.Vivulskio str. 12D, 03221, Vilnius, Lithuania
www.kid-man.com | info@kid-man.com
Tel. +370 5 265 0000

Importator / Distribuitor RO : Kid-Man Aid SRL
Strada Theodor Aman Nr 2-2A, Braila Tel: +40 753 660066

Biodegradable Nitrile Gloves



MAIN FEATURES

- Biodegradable gloves
- Biodegradation within 270 days
- Contains innovative biodegradable technology
- Powder free



Color:
Blue

Type:
Biodegradable

Recommendations for use:



Basic Procedures
and wards



Cleaning



Food Handling

MEDICAL	2017/745	CLASS I	EN 455	AQL 1.5
PPE	2016/425	CATEGORY III	EN 374	AQL 1.5

MAIN INFORMATION

SIZE	REF	PALM WIDTH	OVERALL LENGHT	THICKNESS		
				PALM	CUFF	FINGER
XS	NIGLXS100BDG	76 ± 4 mm	≥240 mm	0.07 ± 0,02 mm	0.06 ± 0,03 mm	0.10 ± 0.03 mm
S	NIGLS100BDG	86 ± 4 mm				
M	NIGLM100BDG	98 ± 4 mm				
L	NIGLL100BDG	107 ± 4 mm				
XL	NIGLXL100BDG	115 ± 4 mm				

MATERIAL	EXTERIOR	INSIDE TREATMENT	FORCE AT BREAK		SHELF TIME	NATURAL LATEX FREE	ACCELERATOR FREE
			Before ageing	After ageing			
Nitrile Latex	Fingertip textured	Chlorinated	≥6N		5 years	Yes	No

Biodegradable Nitrile Gloves

LABELLING

- Each unit displays the following information: product reference, batch number and expiry date.
- Symbols used in compliance with standard EN ISO 15223: Symbols to be used with medical device labels, labelling and information to be supplied.

CE 2777

- Notified Body responsible for certification and ongoing conformity: SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland



EN ISO 374-5:2016
Protects against viruses

- Resistance to bacteria and fungi - pass. Resistance to virus - pass.*The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

PACKAGING / LOGISTIC INFORMATION


AVAILABLE SIZES	PCS/DISPENSER	PCS/TRANSPORT CARTON	DISPENSER	TRANSPORT CARTON
XS	100 pcs.	1000 pcs.	 4 779031 633476	 4 779031 633483
S	100 pcs.	1000 pcs.	 4 779031 633490	 4 779031 633506
M	100 pcs.	1000 pcs.	 4 779031 633513	 4 779031 633520
L	100 pcs.	1000 pcs.	 4 779031 633537	 4 779031 633544
XL	100 pcs.	1000 pcs.	 4 779031 633551	 4 779031 633568



100 PCS PER BOX



10 PCS PER BOX



840 BOXES PER PALLETTE

Sales Unit: Transport Carton

Storage conditions:



UAB Kasko Group | Kid-Man®


A.Vivulskio str. 12D, 03221, Vilnius, Lithuania
www.kid-man.com | info@kid-man.com
Tel. +370 5 265 0000

Importator / Distribuitor RO : Kid-Man Aid SRL
Strada Theodor Aman Nr 2-2A, Braila Tel: +40 753 660066

KID-MAN GRIP gloves are designed for professionals working in demanding environments. Made from flexible and durable nitrile, these gloves are twice as thick as standard thin nitrile gloves, offering enhanced protection. The unique diamond 3D texture ensures a secure grip, even when handling wet or oily items.

MAIN FEATURES

- **Strong grip** – the diamond texture provides stability in greasy and wet environments.
- **Comfortable wear** – inner texture reduces hand sweating.
- **High protection** – shields against select chemicals, bacteria, viruses, and fungi.
- **Durability** – less wear and tear than standard nitrile gloves.
- **Food-safe** – approved for food processing applications.



Color:

Black

Recommendations for use:



Heavy duty



Extra thick



Vehicle repair



Professional service

MEDICAL	2017/745	CLASS I	EN 455	AQL 1.5
PPE	2016/425	CATEGORY III	EN 374	AQL 1.5

MAIN INFORMATION

SIZE	REF	PALM WIDTH	OVERALL LENGTH	THICKNESS		
				PALM	CUFF	FINGER
M	NIGLM50GRIP	95 ±10 mm	≥240 mm	0.17 ± 0,03 mm	0.09 ± 0,03 mm	0.21 ± 0.03 mm
L	NIGLL50GRIP	110 ±10 mm				
XL	NIGLXL50GRIP	120 ±10 mm				
XXL	NIGLXXL50GRIP	130 ±10 mm				

MATERIAL	EXTERIOR	INSIDE TREATMENT	FORCE AT BREAK		SHELF TIME	NATURAL LATEX FREE	ACCELERATOR FREE
			Before ageing	After ageing			
Nitrile	Diamond textured	Chlorinated	≥6N		5 years	Yes	No

LABELLING

- Each unit displays the following information: product reference, batch number and expiry date.
- Symbols used in compliance with standard EN ISO 15223: Symbols to be used with medical device labels, labelling and information to be supplied.

CE 2777

- Notified Body responsible for certification and ongoing conformity: SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland



EN ISO 374-5:2016
Protects against viruses

- Resistance to bacteria and fungi - pass. Resistance to virus - pass.*The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

PACKAGING / LOGISTIC INFORMATION

AVAILABLE SIZES	PCS/DISPENSER	PCS/TRANSPORT CARTON	DISPENSER	TRANSPORT CARTON
M	50 pcs.	500 pcs.	 4 779031 633971	 4 779031 633988
L	50 pcs.	500 pcs.	 4 779031 634008	 4 779031 633995
XL	50 pcs.	500 pcs.	 4 779031 634015	 4 779031 634022
XXL	50 pcs.	500 pcs.	 4 779031 634046	 4 779031 634039



50 PCS PER BOX



10 PCS PER BOX



800 BOXES PER PALLETTE

Sales Unit: Transport Carton

Storage conditions:



UAB Kasko Group | Kid-Man®

A.Vivulskio str. 12D, 03221, Vilnius, Lithuania
www.kid-man.com | info@kid-man.com
Tel. +370 5 265 0000

Importator / Distribuitor RO : Kid-Man Aid SRL
Strada Theodor Aman Nr 2-2A, Braila Tel: +40 753 660066

MAIN FEATURES

- Latex gloves are made of natural rubber (NR) latex, which offers good resistance against many acids and alkalis
- The best fit and feel, latex gloves are very comfortable to wear
- A disadvantage of natural latex is that its natural proteins can cause or trigger latex allergies
- Protein content does not exceed 50 µg/g



Color:
White

Recommendations for use:



Basic Procedures
and wards



Cleaning



Food Handling

MEDICAL	2017/745	CLASS I	EN 455	AQL 1.5
PPE	2016/425	CATEGORY III	EN 374	AQL 1.5

MAIN INFORMATION

SIZE	REF	PALM WIDTH	OVERALL LENGHT	THICKNESS		
				PALM	CUFF	FINGER
XS	LTGLXS100	74 ± 4 mm	≥240 mm	0.10 ± 0,02 mm	0.08 ± 0,03 mm	0.12 ± 0.03 mm
S	LTGLS100	84 ± 4 mm				
M	LTGLM100	94 ± 4 mm				
L	LTGLL100	105 ± 4 mm				
XL	LTGLXL100	114 ± 4 mm				

MATERIAL	EXTERIOR	INSIDE TREATMENT	FORCE AT BREAK		SHELF TIME	NATURAL LATEX FREE	ACCELERATOR FREE
			Before ageing	After ageing			
Natural Rubber Latex	Fingertip textured	Chlorinated	≥6N		5 years	No	No

LABELLING

- Each unit displays the following information: product reference, batch number and expiry date.
- Symbols used in compliance with standard EN ISO 15223: Symbols to be used with medical device labels, labelling and information to be supplied.

CE 2777

- Notified Body responsible for certification and ongoing conformity: SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland



EN ISO 374-5:2016
Protects against viruses

- Resistance to bacteria and fungi - pass. Resistance to virus - pass.*The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

PACKAGING / LOGISTIC INFORMATION

AVAILABLE SIZES	PCS/DISPENSER	PCS/TRANSPORT CARTON	DISPENSER	TRANSPORT CARTON
XS	100 pcs.	1000 pcs.		
S	100 pcs.	1000 pcs.		
M	100 pcs.	1000 pcs.		
L	100 pcs.	1000 pcs.		
XL	100 pcs.	1000 pcs.		



100 PCS PER BOX



10 PCS PER BOX



800 BOXES PER PALLETTE

Sales Unit: Transport Carton

Storage conditions:



UAB Kasko Group | Kid-Man®

A.Vivulskio str. 12D, 03221, Vilnius, Lithuania
www.kid-man.com | info@kid-man.com
Tel. +370 5 265 0000

Importator / Distribuitor RO : Kid-Man Aid SRL
Strada Theodor Aman Nr 2-2A, Braila Tel: +40 753 660066

MAIN FEATURES

- Thin powder-free vinyl glove, ambidextrous and with rolled rim for extra cuff strength
- Good viral barrier properties, protects against variety of chemicals
- Easy to don and remove
- Free of latex, powder and proteins, good alternative to those suffering from type I allergies
- The gloves are suitable for all food contact



Color:
Clear

Recommendations for use:



Basic Procedures
and wards



Cleaning



Food Handling

MEDICAL	2017/745	CLASS I	EN 455	AQL 1.5
PPE	2016/425	CATEGORY III	EN 374	AQL 1.5

MAIN INFORMATION

SIZE	REF	PALM WIDTH	OVERALL LENGTH	THICKNESS		
				PALM	CUFF	FINGER
S	VIGLS100	85 ± 5 mm	≥240 mm	Min 0.08 mm	Min 0.06 mm	Min 0.10 mm
M	VIGLM100	95 ± 5 mm				
L	VIGLL100	105 ± 5 mm				
XL	VIGLXL100	115 ± 5 mm				

MATERIAL	EXTERIOR	INSIDE TREATMENT	FORCE AT BREAK		SHELF TIME	NATURAL LATEX FREE	ACCELERATOR FREE
			Before ageing	After ageing			
Polyvinyl chloride (PVC)	Smooth	Powder Free	Min. 3.6N		5 years	Yes	No

LABELLING

- Each unit displays the following information: product reference, batch number and expiry date.
- Symbols used in compliance with standard EN ISO 15223: Symbols to be used with medical device labels, labelling and information to be supplied.

CE 2777

- Notified Body responsible for certification and ongoing conformity: SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland



EN ISO 374-5:2016
Protects against viruses

- Resistance to bacteria and fungi - pass. Resistance to virus - pass.*The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

PACKAGING / LOGISTIC INFORMATION

AVAILABLE SIZES	PCS/DISPENSER	PCS/TRANSPORT CARTON	DISPENSER	TRANSPORT CARTON
S	100 pcs.	1000 pcs.	 4 779031 633148	 4 779031 633193
M	100 pcs.	1000 pcs.	 4 779031 633155	 4 779031 633209
L	100 pcs.	1000 pcs.	 4 779031 633162	 4 779031 633216
XL	100 pcs.	1000 pcs.	 4 779031 633179	 4 779031 633223



100 PCS PER BOX



10 PCS PER BOX



800 BOXES PER PALLETTE

Sales Unit: Transport Carton

Storage conditions:



UAB Kasko Group | Kid-Man®

A.Vivulskio str. 12D, 03221, Vilnius, Lithuania
www.kid-man.com | info@kid-man.com
Tel. +370 5 265 0000

Importator / Distribuitor RO : Kid-Man Aid SRL
Strada Theodor Aman Nr 2-2A, Braila Tel: +40 753 660066

MAIN FEATURES

- Light nitrile construction: thin, elastic and flexible for a comfortable fit without reducing sensitivity
- Powder-free, non-sterile, ambidextrous with beaded cuff for added cuff strength
- Textured fingertips support secure wet and dry grip while maintaining excellent tactile sensitivity
- High-level barrier protection against many chemicals and infectious agents, especially viruses
- Compatible with water and alcohol based solutions and gels
- The gloves are suitable for all food contact



Color:
Blue

Recommendations for use:



Basic Procedures
and wards



Cleaning



Food Handling

MEDICAL	2017/745	CLASS I	EN 455	AQL 1.5
PPE	2016/425	CATEGORY III	EN 374	AQL 1.5

MAIN INFORMATION

SIZE	REF	PALM WIDTH	OVERALL LENGTH	THICKNESS		
				PALM	CUFF	FINGER
XS	NIGLXS100LIGHT	76 ± 4 mm	≥240 mm	0.07 ± 0.02 mm	Min 0.07 mm	0.09 ± 0.02 mm
S	NIGLXS100LIGHT	86 ± 4 mm				
M	NIGLXS100LIGHT	98 ± 4 mm				
L	NIGLXS100LIGHT	107 ± 4 mm				
XL	NIGLXS100LIGHT	115 ± 4 mm				

MATERIAL	EXTERIOR	INSIDE TREATMENT	FORCE AT BREAK		SHELF TIME	NATURAL LATEX FREE	ACCELERATOR FREE
			Before ageing	After ageing			
Nitrile	Fingertip	Chlorinated, powder free	≥6N		5 years	Yes	No

LABELLING

- Each unit displays the following information: product reference, batch number and expiry date.
- Symbols used in compliance with standard EN ISO 15223: Symbols to be used with medical device labels, labelling and information to be supplied.

CE 2777

- Notified Body responsible for certification and ongoing conformity: SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland



EN ISO 374-5:2016
Protects against viruses

- Resistance to bacteria and fungi - pass. Resistance to virus - pass.*The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

PACKAGING / LOGISTIC INFORMATION

AVAILABLE SIZES	PCS/DISPENSER	PCS/TRANSPORT CARTON	DISPENSER	TRANSPORT CARTON
XS	100 pcs.	1000 pcs.	 4 779031 635036	 4 779031 635043
S	100 pcs.	1000 pcs.	 4 779031 635067	 4 779031 635050
M	100 pcs.	1000 pcs.	 4 779031 635074	 4 779031 635081
L	100 pcs.	1000 pcs.	 4 779031 635098	 4 779031 635104
XL	100 pcs.	1000 pcs.	 4 779031 635111	 4 779031 635128



100 PCS PER BOX



10 PCS PER BOX



1080 BOXES PER PALLETTE

Sales Unit: Transport Carton

Storage conditions:



UAB Kasko Group | Kid-Man®

A.Vivulskio str. 12D, 03221, Vilnius, Lithuania
www.kid-man.com | info@kid-man.com
Tel. +370 5 265 0000

Importator / Distribuitor RO : Kid-Man Aid SRL
Strada Theodor Aman Nr 2-2A, Braila Tel: +40 753 660066

Regulatory & requirements

What do the pictograms and symbols on the glove boxes mean?

	Manufacturer		Date of manufacture		Temperature limitation		Consult instructions for use
	Lot / batch number		Expiry date		For single use only		Recycling packaging
	Product number		Keep dry		Non-sterile		This side up
	Medical device		Keep away from sunlight		Suitable for food contact		Handle with care
	Serial number						

Let's get in touch

 UAB Kasko Group | Kid-Man®
 A.Vivulskio str. 12D
 03221, Vilnius, Lithuania
 Tel. +370 5 265 0000
order@kid-man.com

 Kid-Man AID S.R.L
 Str. Theodor Aman nr.2-2A
 810031 Braila, Romania
 Telefon +40 753 660066
ro.order@kid-man.com



www.kid-man.com