

STEELMAN Start | 04-020

MAIN FEATURES

- Universal wheelchair with strong and reliable construction is intended for people weighing up to 120 kg with temporary or long-term mobility problems.
- The wheelchair can be used both indoors and outdoors, to move independently by rotating metal driving arcs located on the rear wheels.
- The wheelchair can be driven by the attendant using pushing handles mounted above the backrest.
- A number written on the round sticker on the armrest (e.g. 46) indicates that the width of the wheelchair seat is 46 cm.
- The structure with a cross-shaped element, so the wheelchair is easily folded it is convenient to transport it.
- Width of the folded wheelchair – only 26 cm, the height of the backrest – 40.5 cm, the depth of the seat – 41 cm.
- The seat and backrest of the wheelchair are padded.
- Wheelchair have parking brakes.
- Armrests are folded and removed using fasteners, the distance between the armrests and the seat is 22 cm.
- Removable, rotating to the sides, height-adjustable footrests with removable calf support strap.
- Fast installation rear wheels.
- The wheelchair comes with reflectors, adjustment key, instructions for use and maintenance in Lithuanian.



Color:
Black



Recommendations for use:



For nursing



For indoors and outdoors

Care:

It is recommended to store the wheelchair in a dry place protected from direct sunlight, to wet clean it periodically without using aggressive cleaners, or to cover it.

Warranty:

The device has a 24-month warranty (excluding front and rear wheels) from the date of purchase.

MAIN INFORMATION

TYPE	MORE SIZES	WEIGHT, KG	WHEEL MATERIAL	PRODUCT MATERIAL	ACCESSORIES /FURNITURE-MATERIAL
04-020-2	36	17.6	PU front casters; Pneumatic rear wheel.	The structure is made of steel, painted; Driving arcs and rims – aluminium; PU armrests with alliuminum core inside;	On the edge of the pocket, legrest strap-velcro straps.
	38	17.8			
	41	18			
	43	18.2			
	46	18.4			
04-020-3	48	18.6	PU front casters; Solid rear wheel.	Nylon seat- strong PVC inside; Rubber plastic handles.	
	38	18			
	41	18.4			
	43	18.8			
	46	19.2			
	48	19.6			

MAX LOAD, KG	FOOTREST POSITIONS	GENERAL HEIGHT, CM	GENERAL WIDTH, CM	GENERAL LENGTH, CM	BACK WHEEL SIZE,CM	FRONT WHEEL SIZE, CM
120	6	92	Seat width + 20	107	61	20 x 5

PACKAGING LOGISTIC INFORMATION

Package info:

REF	SIZE	EAN-13	PCS PER BOX	PCS PER PALETTE
04-020-2	36		1 pc	8 pcs
	38			
	41			
	43			
	46			
	48			

04-020-3	38		1 pc	8 pcs
	41			
	43			
	46			
	48			

Package info:

REF	04-020-2 / 04-020-3	QTY	NW	GW	MAX LOAD	SEAT WIDTH
DIM	L800 x W290 x H930 mm	1 pc	18 kg	20 kg	120 kg	Label with seat size attached on the box and the wheelchair frame

Storage conditions:



The quality management system of Kasko Group
UAB meets the standards ISO 9001 and ISO 14001

Manufacturer:

UAB Kasko Group | Kid-Man®
A.Vivulskio str. 12D, 03221, Vilnius, Lithuania
www.kid-man.com | info@kid-man.com
Tel. +370 5 265 0000

Importer / Distributor RO : Kid-Man Aid SRL
Strada Theodor Aman Nr 2-2A, Braila Tel:Tel: +40 753 660066